



V1.0

Test instructions

OpenCASA




Table of Contents


Introduction	1
Simulation	1
Chemotaxis	2
Dataset	2
Analyze file	2
Simulated experiment	3
Ch-Index	4
Bootstrapping	5
Experiment analysis	6
Motility	7
Analyze file	7
Analyze directory	8
Analyze multiple directories	8
Viability	9
Morphometry	11
Validation data	12
Motility	12
Viability	13
Morphometry	14

Introduction

In this document are presented the instructions to test the software. Except for the simulation module, all configuration parameters of OpenCASA that were set to carry out each analysis, are available at the corresponding directory of the test dataset that goes with this document.

Simulation

This module generates a simulation of 100 cells in movement. There are three parameters needed to generate a simulation:




Beta equals 0 means the population will behave without any bias in the directionality of the swimming path. Values greater/lower than 0 will force the cells to move to the right/left respectively. Responsiveness is the percentage of cells that will be forced to some direction when beta is not zero. The length of the simulation is the duration (in frames) of the simulated video. In order to test visually the effect and aspect of the simulations, check the following situations:

Condition	Parameters	Expected response
Control	Beta=0; Responsiveness = indifferent	Random swimming paths in all directions
Attraction to the right (moderate)	Beta = 2 Responsiveness = 50	The half of the population will eventually move to the right side of the screen.
Attraction to the right (High level)	Beta = 2 Responsiveness = 100	All the population will eventually move to the right side of the screen.
Attraction to the left (moderate)	Beta = -2 Responsiveness = 50	The half of the population will eventually move to the left side of the screen.

Chemotaxis

The module allows the user to carry out three types of analysis: analyze a file, an experiment, or simulate an experiment. Each analysis has different purposes: analyze a single file allows to check the distributions of the instantaneous displacements of the whole population, and also to take a look at the relative trajectories; Experiment analysis tries to characterized the bias in the swimming direction of a particular set of sperm populations, comparing it to a control set; and Simulated Experiment helps the user to validate and learn the basis of an Experiment analysis.



Dataset

In order to test the module, various videos are provided. In the chemotaxis folder of the validation data, it is possible to find the "Experiment" folder, with 4 simulations: two about control conditions, and two of biased populations.

Analyze file

Click on *Chemotaxis Module* => *Analyze File*. After selecting a video, the program will show two pictures: one about the relative trajectories detected by the module, and the other about the distribution of the instantaneous directionality angles of the whole population. In color green is drawn the chemotactic cone, and in color red the chemotaxis direction set by the user.




Figure 1 - Relative trajectories detected by the software. Each trajectory is drawn as if it started from the origin of coordinates.





Figure 2- Differences on the distributions of the instantaneous directionality angles.
(a) a control sample; (b) a biased sample

Simulated experiment


This analysis emerge from the need of validate the bootstrapping analysis for a large set of files. The way the program analyze a simulated set of videos is the same that is carried out with a real ones, but the process of creating the simulations is embedded in order to save time to the user. Other way would be to create manually a set of videos using the simulation module, place them in a folder (correctly structured as it is explained in the user's manual), and later carry out an "Experiment analysis".

To simulate an experiment: click on "**Chemotaxis --> Simulated Experiment**". There will appear two options: Ch-Index or Bootstrapping.

Test instructions



Each option is explained below but for both cases, the user will have to set the parameters of the simulations. As in the Simulation module, Beta equals 0 means the population will behave without any bias in the directionality of the swimming path. Values greater or lower than 0 will force the cells to move to the right or left respectively. Responsiveness is the percentage of cells that will be forced when beta is not zero. The length of all simulations in this module is fixed to 500 frames. The number of simulations means the number of pairs control-biased samples that will be simulated. For example, a value of 10 will create 10 simulations with Beta equals 0, paired up with another 10 simulations created with those parameters set by the user (in total, 20 simulations).



Ch-Index

The Ch-Index analysis gives a measure of the percentage of the instantaneous directionality angles that point to the gradient direction respect to those not pointing to that gradient (an explanation in more detail can be found in the attached research article). In order to validate this analysis, it is necessary to check that, in control conditions (no bias in the directionality), the percentage of the instantaneous displacements pointing to any direction has to be proportional to the area taken into account. For example, if the user defines a chemotactic cone of 60 degrees, if there is no bias in the directionality, statistically the Ch-Index should be $\frac{60^\circ}{360^\circ} * 100 \approx 16,67\%$. Else, if there exists a bias in the chemotaxis direction, the percentage should be greater. This test can be seen in Figure 3, setting Beta equals 2 and Responsiveness equals 50%.

Test instructions






Figure 3 - Ch-Index analysis of 4 pairs of simulations. The biased ones were created with this parameters: Beta=2; Responsiveness=50%. The Ch-Index was calculated taking into account all directions on the instantaneous displacements (option 1, see paper).

However, there is an option to take into account only the displacements in the same and opposite direction of the chemotaxis gradient. In that case, if the user defines a chemotactic cone of 60 degrees, if there is no bias in the directionality, statistically the Ch-Index should be $\frac{60^\circ}{60^\circ + 60^\circ} * 100 = 50\%$. This option can be set in the chemotaxis tab of the settings menu:



Analysing a set of simulations with the same parameters as the previous example, but comparing in this case only opposite directions, we get the following results:



Bootstrapping

The Bootstrapping analysis act as a binary classifier (positive/negative) looking at the distribution of the instantaneous displacements of two populations: a control condition and a hypothesized biased one. The software will label a sample as positive if, statistically compared to a control condition, exists a bias in the direction of the chemotaxis gradient. By definition, this clasification has a 5% rate of false positives, so obtaining percentage of positives upper this value would indicate a chemotaxis phenomena in our samples. To test this analysis, we can carry out two different analysis. First, we can test if, analysing a set of NON-biased samples, we obtain only a 5% of false positives (take into account that this


Test instructions

percentage in the theoretical one obtained when the sample size is infinite, so a little variance around this value could be expected). To do it, we can simulate for example, 100 samples with Beta equals 0, and count the number of positives that we obtain. It is recommended to do it several times in order to be sure that this number is close to the theoretical 5%. Later, it is necessary to test that this percentage grows up when a bias is introduced. To do that, just simulate again setting the value of Beta above 0 (for example, Beta = 0.5). The value of Responsiveness is also important because is the number of cells that actually respond to the chemotactic effect. For the test below, Responsiveness has been set to 50%.

Results from Simulation (Beta: 0.5, Responsiveness: 0.5)						
File Edit Font						
ID	OR	Threshold	Result	Type	Source	COD
0	1.344	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\0	NO
1	1.366	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\1	NO
2	1.362	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\2	NO
3	1.293	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\3	NO
4	1.196	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\4	NO
5	1.395	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\5	NO
6	1.450	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\6	NO
7	1.527	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\7	NO
8	1.334	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\8	NO
9	1.374	1.116	POSITIVE	beta: 0.5; resp: 0.5	Beta: 0.5; Resp: 0.5\9	NO


Experiment analysis

In order to test the third analysis, four videos are provided. After clicking on the button “Experiment Analysis” a dialog window is opened. After choosing between Ch-Index or Bootstrapping analysis, select the folder “Experiment” inside the Chemotaxis directory and the analysis will be carried out as in the previous section for simulated data.




Motility

In order to test the motility module, a set of videos is provided. This module allows the user to analyze files at three levels: one single file, a directory with different videos, or a directory with multiple directorie.



Analyze file

When a single file is analyzed, the original video is shown with the detected trajectories overwriting it. The user can check the correspondences between the trajectories and the cells, using the slider at the bottom of the window:



Slider

Also, the analysis return two reports: one for the individual analysis, and other with the average values for the whole sample.

Test instructions

File	Edit	Font	Track Number	Length (frames)	VSL (um/s)	VCL (um/s)	VAP (um/s)
			1	55	9.715	43.776	9.439
			2	55	88.519	280.609	108.99
			3	50	28.374	111.109	40.886
			5	16	53.286	420.600	156.08
			9	50	6.696	178.982	33.249
			13	59	28.708	189.556	105.29
			14	34	58.891	206.767	96.448
			17	17	178.218	294.511	187.00
			19	59	69.680	222.167	90.536
			21	59	32.152	96.709	41.413

File	Edit	Font	Motile trajectories	VSL Mean (um/s)	VCL Mean (um/s)	VAP Mean (um/s)
			99	103.441	237.052	125.623

Analyze directory

The directory analysis is just the automation of the file analysis but for multiple files. The only difference is that in this analysis, the trajectories over the videos are not shown.


Analyze multiple directories

This analysis is the automation for the previous analysis. Selecting this option, multiple directories are analyzed, but in this case, the report given by the program shows only the average values calculated for each subdirectory, independently of how many samples has each directory. The report shows the average values (one directory for each row):

File	Edit	Font	Motile trajectories	VSL Mean (um/s)	VCL Mean (um/s)	VAP Mean (um/s)	LIN Mean	WOB Mean	STR Mean	ALH_Mean Mean (um)	ALH_Max Mean (um)
			878	101.566	222.102	120.770	45.591	54.060	83.101	0.965	2.711
			694	129.891	258.460	144.940	49.877	55.990	88.456	0.951	2.621
			847	139.060	270.739	155.535	51.923	57.923	88.975	0.845	2.574
			764	131.323	283.169	150.029	46.341	52.738	87.461	1.032	2.983
			559	144.688	280.689	160.997	52.757	58.217	90.564	0.849	2.498
			647	126.039	301.753	151.398	42.262	50.335	83.753	1.301	3.465
			838	129.749	294.371	154.585	44.400	52.629	83.621	1.161	3.357
			714	125.415	285.017	150.239	44.651	52.852	83.977	1.109	3.154
			952	123.507	308.789	154.984	40.578	50.746	79.299	1.368	3.618
			548	113.394	268.855	134.466	42.286	50.101	83.480	1.315	3.136
			944	128.953	265.337	153.214	49.504	58.436	84.171	0.882	2.751
			769	145.902	298.288	168.814	50.232	57.700	86.601	1.008	2.920
			624	116.481	237.155	132.625	50.517	57.056	87.968	0.906	2.447
			622	154.042	285.610	171.260	55.084	60.756	90.391	0.861	2.556
			984	120.914	261.802	143.310	47.299	55.795	83.850	0.971	2.907
			748	153.954	301.295	175.085	51.810	58.493	88.161	0.884	2.786
			556	120.589	281.343	141.454	43.698	50.729	85.314	1.237	3.188
			446	102.268	252.139	129.420	40.006	51.666	76.623	1.258	3.156
			541	134.759	276.874	154.866	48.812	55.941	86.548	1.113	3.026
			352	132.065	273.905	148.568	49.120	55.315	87.834	0.951	2.678
			803	125.805	258.462	151.457	48.583	58.616	81.542	0.937	2.743
			291	137.605	246.392	151.353	56.607	62.109	91.031	0.721	2.139
			292	119.392	250.557	139.429	47.541	55.716	84.652	1.076	2.856
			385	151.285	259.436	165.634	57.959	63.556	90.825	0.734	2.254
			693	121.535	258.883	141.386	46.950	54.772	85.182	1.087	2.979
			337	135.153	264.210	146.159	52.293	56.541	92.303	0.872	2.370
			564	113.009	230.760	126.576	50.660	56.484	88.896	0.790	2.209
			819	136.522	298.158	162.897	46.046	54.656	83.768	1.084	3.180
			495	139.248	283.827	157.752	49.723	56.181	88.010	0.978	2.884
			742	123.512	305.011	149.388	40.891	49.268	82.041	1.638	3.802
			771	86.981	214.308	108.256	40.808	50.545	80.610	1.301	3.099
			568	97.834	274.217	128.243	37.089	48.260	75.611	1.518	3.718
			457	112.153	277.145	134.738	40.581	48.938	81.456	1.401	3.366


Viability

In order to test the viability module, a set of images is provided. This module allows the user to choose between a single file analysis (only one image), or a directory analysis (all images inside the selected directory). Both analysis are equivalent so here we will explain only how to carry out a directory analysis.




Once the directory has been selected, the program loads the first image and it tries to classify the cells inside the image between viable and non-viable cells. In case that some cells are labeled in a wrong way, the user can modify the thresholds used for red and green channels, just moving the corresponding slider:

Sliders



Test instructions

Once the user considers that the classification is correct, clicking on bottom-right button “Save and Next Image”, the results will be added to the report window. In this window, on the left we will see the results for each analyzed image, and on the right, the absolute and relative number of analyzed cells. This second part is updated each time a new image is analyzed.




The screenshot shows a window titled "Viability results" with a menu bar (File, Edit, Font) and a table of data. The table has columns for Total, Viable, Non-Viable, % Viable, % Non-Viable, Sample, and Filename. The data is organized into rows for different samples, with the first row showing a total of 326 cells, 210 viable, and 116 non-viable.


Total	Viable	Non-Viable	% Viable	% Non-Viable	Sample	Filename	Total Cells	Total Viable	Total Non-Viable	% Total Viable	% Total Non-Viable
71	54	17	76.056	23.944	Muestra 12	12 (2)_01	326	210	116	64.417	35.583
90	52	38	57.778	42.222	Muestra 12	12 (2)_02					
79	51	28	64.557	35.443	Muestra 12	12 (2)_03					
86	53	33	61.628	38.372	Muestra 12	12 (2)_04					

Morphometry

In order to test the morphometry module, a set of images is provided. This module allows the user to choose between a single file analysis (only one image), or a directory analysis (all images inside the directory). Both analysis are equivalent so here we will explain how to carry out a directory analysis.



Once the directory has been selected, the program loads the first image and it tries to determine the optimal threshold to detect the contour of the cell. In this module, the user has to select manually those cells that have been well outlined. When it occurs, the morphometric values of the cell are added to the report window. If a cell has not been well outlined, the user can adjust manually the threshold just moving the slider on the top of the window.



Validation data

Considering that the simulation module can be tested visually, the chemotaxis module is perfectly validated with these simulations, and the rest of the modules have been validated using real data, the results of that validation are provided in order to allow the user to replicate and confirm the analysis. For that reason, a subset of the validation data is provided together with the results of the whole analysis. Taking into account the type of analysis, and the identification of each sample, the user has to be able to replicate the corresponding results.

Motility

The file “Motility results.xlsx” contains the results of the analysis carried out for each sample (both for the ISAS and OpenCASA software). Also, two samples with 8 videos each sample are provided. The steps that the user has to follow to replicate the results of the analysis of that samples are:

1. Set the OpenCASA parameters as said in the file “OpenCASA Parameters Motility.txt” placed in the motility folder.
2. Select the module Motility and click on “Multiple directories”
3. Select the folder “data” contained in the motility folder

After that, compare that the results obtained in the analysis are the same that those that figure on the excel file. To identify each sample, just move to the end of the report to check the last part of the filename and look for the sample ID on the first column of the spreadsheet:

Sample	Motile trajectories	VSL Mean (um/s)	VCL Mean (um/s)	VAP Mean (um/s)	UN Mean	WOB Mean	STR Mean	ALH Mean Mean (um)	ALH Max Mean (um)	BCF Mean (Hz)	DAI
1	878	101.566	222.102	120.77	45.591	54.06	83.101	0.965	2.711	36.514	
2	694	129.891	258.46	144.94	49.877	55.99	88.456	0.951	2.621	41.203	
3									2.574	45.201	
4									2.983	42.143	
5									2.498	43.224	
6									3.465	40.355	
7									3.357	40.932	
8									3.154	40.398	
9									3.618	39.095	
10									3.136	36.611	
11									2.751	42.337	
12									2.920	43.377	
13									2.447	39.91	
14									2.556	43.226	
15									2.907	41.822	
16									2.785	44.726	
17									3.188	37.836	
18									3.156	34.558	
19									3.026	40.215	
20									2.678	39.069	
21									2.743	39.401	
22									0.721	44.016	
23									1.076	38.059	
24									0.734	43.985	
25									1.087	38.692	
26									0.872	41.478	
27									0.790	40.944	
28									1.084	40.883	
29									0.978	41.699	
30									1.638	36.118	
31									1.301	31.119	
32									1.518	31.274	
33									1.401	35.258	
34									1.112	35.152	
35									2.510	35.338	

Viability

The file “Viability results.xlsx” contains the results of the analysis carried out for each sample. Also, five samples with 52 images in total are provided. The steps that the user has to follow to replicate the results of the analysis of that samples are:

1. Set the OpenCASA parameters as said in the file “OpenCASA Parameters Viability.txt” placed in the viability folder.
2. Select the module Viability and click on “analyze directory”
3. Select the folder of one sample contained in the data folder.
4. Analyze the images of that sample and save the results
5. Repeat that process for each sample directory
6. After analyzing all samples, compare the values obtained with those of the same samples provided in the Viability results.xlsx file

The validation results

The screenshot displays an Excel spreadsheet titled 'Resultados muestras 5-36'. The main table lists 36 samples with columns for Sample, % viables, % non-viables, and Total cells. A red arrow points to the 'Total cells' column. To the right, a summary table for 'Muestra 7' compares user results with the file's results.

Total cells	
mean	std
670.551724	308.493411


	average of %	% with absolute values	Total
Viabiles	67.9627	67.7852349	505
Non-Viables	32.0373	32.2147651	240
			745

Morphometry

The file “Morphometry results.xlsx” contains the results of the analysis carried out for each sample (both for the ISAS and OpenCASA software). Also, three samples with 100 images in total are provided. The steps that the user has to follow to replicate the results of the analysis of that samples are:

1. Set the OpenCASA parameters as said in the file “OpenCASA Parameters Morphometry.txt” placed in the viability folder.
2. Select the module Morphometry and click on “analyze directory”
3. Select the folder of one sample contained in the data folder.
4. Analyze the images of that sample and save the results. Select around 200 cells/sample, equally distributed across all images.
5. Repeat that process for each sample directory
6. After analyzing all samples, compare the values obtained with those of the same samples provided in the Morphometry results.xlsx file.

The validation results



Sample ID	Total Cells	Area	Mean gray value	PERIM	FERET (L)	Min_Feret (W)
001-97690	202	30.52458911	112.0969851	22.63931683	8.287950495	4.626
002-81302	204	31.76346569	154.0151618	22.90215196	8.520553922	4.698
003-16384	204	30.96740686	156.232652	22.14803431	8.052196078	4.83
004-7170	204	30.65879902	126.9159755	22.69284804	8.448622549	4.598
005-97689	204	30.88681373	119.8051422	22.60832353	8.277617647	4.714
006-16375	204	31.93037255	111.265549	22.90313235	8.333818627	4.816
007-61871	204	31.9500049	139.5103431	23.01137745	8.498985294	4.687
008-16383	204	30.8322402	99.83638235	22.65676961	8.269102941	4.722
009-5347	205	31.13039512	128.1255707	22.58282927	8.31382439	4.773
010-753	200	29.201655	88.61499	21.9796	8.09697	4.5
011-5008	200	29.425945	110.025495	21.844905	8.054535	4.6
012-6085	208	30.27119231	93.92761058	22.35507692	8.267394231	4.633
013-6091	200	30.526355	104.0199	22.50462	8.402215	4.5
014-6125	205	31.55706341	93.25346341	22.84330732	8.414512195	4.742
015-7150	200	29.414505	87.20177	21.984995	8.127985	4.5
016-7151	202	31.35138119	91.36432673	22.71045545	8.328450495	4.761
017-7170	208	30.79721154	87.59507212	22.63678846	8.458913462	4.616
018-5008	204	28.81491667	122.6876765	21.57395098	7.970142157	4.564
020-5348	204	31.38675	111.287049	22.46971078	8.209176471	4.81
021-5350	200	31.268485	119.306	22.53635	8.276185	4.7
022-7161	204	31.071	135.4472353	22.51659314	8.291691176	4.748
023-9201	210	31.17112857	140.3046762	22.59303333	8.364585714	4.694
024-5347	208	30.71009615	102.6293077	22.44624519	8.252206731	4.741
025-5348	208	31.80598558	108.6089471	22.55971635	8.258072115	4.85
026-9201	222	29.66798649	101.7584505	22.07135586	8.18940991	4.563828829
027-97689	216	31.48284722	109.5584213	22.64615278	8.359361111	4.749532407
028-5008	214	29.06092991	121.9877009	21.6689486	8.003242991	4.58235514
029-5348	231	31.05004329	110.3488658	22.36737229	8.173909091	4.787735931
030-5350	210	30.09894286	120.5549714	22.15199048	8.161861905	4.643342857
031-7150	213	29.49221596	91.18520188	21.99075587	8.088267606	4.643685446
032-7151	201	31.25090547	96.65122388	22.53126368	8.278985075	4.784850746
033-9201	210	31.30850952	112.825019	22.60827619	8.362352381	4.714709524
034-80905	211	31.23971564	100.0855261	22.34725592	8.128838863	4.858075829
035-83992	209	31.37808612	104.1360957	22.740311	8.390253589	4.737708134
036-5349	216	30.75215278	98.29307407	22.24647222	8.1765	4.733416667

The results obtained by the user for the same sample